

Clinical Accessories

Trusted. Reliable. Simple.



Contents

PATIENT MONITORING	4	Patient Breathing Circuits	47
ECG	6	Circuit Startup Kits	47
Recording Paper	6	Circuit Kits	48
Trunk Cable Accessories	7	Circuits	49-50
Trunk Cables	8-9	Circuit Components	51
Leadwires	10-13	Flow Sensors – Respiratory	52
Electrodes	14-15	Respiratory Flow Sensors and	
Temperature	16	Expiratory Valve Assemblies	52
Cables	16	Cables	53
Probes	17	Heater	53
Invasive Blood Pressure	18	Flow Sensors – Anesthesia	53
Cables	18-19	Nebulizer Accessories	54
Cardiac Output	20	Anesthesia O2 Cells	54
Cables	20, 23	Respiratory Safety Guard	55
Probes	21-22	Pressure Catheters	55
Other Accessories	22	DIAGNOSTIC CARDIOLOGY	56
Catheter Kits	23	MAC™, CASE™ and CardioSoft™	58
Pulse Oximetry	24	ECG	58
Sensors	24-31	Cables	58-61
Cables	25, 27, 30	Leadwires	58-61
Tape & Wrap	26, 28, 31	Accessories	61
Other Accessories	26	Recording Paper	62
Probes	29	Electrodes	63
Clip	31	SEER® Holter Monitors	64
Adequacy of Anesthesia	32	ECG – CM 4000	64
Cable	32-33	Leadwire Sets	64
Sensors	32-33	Accessories	65
Electrodes	33	ECG – SEER® 12	65
Neurology	34	Leadwire Sets	65
EEG Leadsets	34	Accessories	65
EEG Electrodes	34-35	ECG – SEER® 1000	66
AEP Accessories	35	Leadwire Sets	66
Patient Spirometry	36	Accessories	67
Sensors	36	ECG – SEER® Light	68
Tubes	37	Leadwire and Cable Kits	68
Kits	37	Accessories	69
Gas Monitoring	38		
Airway Adapter Sets	38		
Gas Sampling Lines	39		
Airway Adapters	39		
Dead Space Adapters	40		
Water Traps	40		
Gas Scavenging	41		
Calibration Gases	41		
Regulators	41		
ANESTHESIA & RESPIRATORY CARE	42		
Absorbents	44		
Absorbers	44		
Canisters	44		
Dust Filter	44		
Vaporizer Bottle Adapters	44		
Filters & Breathing Masks	45		
Reservoir Bags & Gas Scavenging	46		
Anesthesia Breathing Bags	46		
Scavenging Bags, Tubes, & Connectors	46		

Welcome

GE Healthcare offers you a simple, reliable one-stop solution for high-quality GE Healthcare-compatible clinical accessories. We provide an expanded inventory of products that have been verified by our engineers. The compatible devices listed may not be available for sale and/or cleared/approved in all countries. Check with your representative for product availability.

PATIENT MONITORING



PATIENT MONITORING

ECG

- Electrodes
- Leadwires
- Recording Paper
- Trunk Cable Accessories
- Trunk Cables

Temperature

- Cables
- Probes

Invasive Blood Pressure

- Cables

Cardiac Output

- Cables
- Catheter Kits
- Other Accessories
- Probes

Pulse Oximetry

- Cables
- Clip
- Other Accessories
- Probes
- Sensors
- Tape & Wrap

Adequacy of Anesthesia

- Cable
- Electrodes
- Sensors

Neurology

- AEP Accessories
- EEG Electrodes
- EEG Leadsets

Patient Spirometry

- Kits
- Sensors
- Tubes

Gas Monitoring

- Airway Adapter Sets
- Airway Adapters
- Calibration Gases
- Dead Space Adapters
- Gas Sampling Lines
- Gas Scavenging
- Regulators
- Water Traps

ECG



ECG RECORDING PAPER



MS ECG Recording Thermal Paper

Length	CARESCAPE™ B450/B650 Monitors, CARESCAPE™ B850 Monitor with PRN 50-M+ printer	B20(i)/B40(i) Monitors	B1x5 Monitors	CARESCAPE™ V100 and VC150 Monitors	Part Number	Quantity
57.2 mm x 24.4 m / 2.25 in x 80 ft				X	2106825-001	10
50 mm x 25.6 m / 2 in x 84 ft	X	X	X		2106824-001	20
50 mm x 27.4 m / 2 in x 90 ft	X	X	X		2106823-001	48

For ECG paper for Diagnostic cardiology and Maternal infant care equipment, please go to pages 62 and 72.

Not available in all countries



2106403-003



2106403-001



2106399-001



2106399-002

ECG Trunk Cable Accessories

Description		Part Number	Quantity
DIN Converter	3-Lead	2106403-003	1
	5-Lead	2106403-001	1
Leadwire separator		2106399-001	3
Bedsheet clip		2106399-002	1

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module E-PSM
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ Monitor with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry¹

Specifications available on request

(1) Not available in the US

Not available in all countries

ECG TRUNK CABLES



ECG Trunk Cables – AHA			
Description	Length	Part Number	Quantity
Neonatal DIN 3-lead	3.6 m/12 ft	2106306-001	1
	1.2 m/4 ft	2106306-002	1
3/5-lead	3.6 m/12 ft	2106305-001	1
	1.2 m/4 ft	2106305-002	1
	6 m/20 ft	2106311-001	1
3/5-lead w/ESU filter	3.6 m/12 ft	2106310-001	1
6-lead	3.6 m/12 ft	2106307-001	1
	1.2 m/4 ft	2106307-002	1
12-lead	3.6 m/12 ft	2106308-001	1
	1.2 m/4 ft	2106308-002	1

Not available in all countries



ECG Trunk Cables – IEC			
Description	Length	Part Number	Quantity
Neonatal DIN 3-lead	3.6 m/12 ft	2106306-003	1
	1.2 m/4 ft	2106306-004	1
3/5-lead	3.6 m/12 ft	2106305-003	1
	1.2 m/4 ft	2106305-004	1
	6 m/20 ft	2106311-002	1
3/5-lead w/ESU filter	3.6 m/12 ft	2106310-002	1
6-lead	3.6 m/12 ft	2106307-003	1
	1.2 m/4 ft	2106307-004	1
12-lead	3.6 m/12 ft	2106308-003	1
	1.2 m/4 ft	2106308-004	1

Not available in all countries

ECG LEADWIRES



2106385-002



2106381-001



2106383-002



2106394-001



2106395-002



2106395-003

ECG Leadwires Set Snap – AHA

Description	Connection	Length	Part Number	Quantity
3-lead	Molded	74 cm/29 in	2106385-001	1
		130 cm/51 in	2106385-002	
5-lead	Non-molded	74 cm/29 in	2106381-001	1
		130 cm/51 in	2106381-002	
		mix 74 cm/29 in, 130 cm/51 in	2106381-005	
	Molded	74 cm/29 in	2106383-001	
		130 cm/51 in	2106383-002	
		mix 74 cm/29 in, 130 cm/51 in	2106383-005	
5-lead V2-6	Non-molded	130 cm/51 in	2106394-001	1
		74 cm/29 in	2106394-003	
6-lead	Non-molded	74 cm/29 in	2106395-001	1
		130 cm/51 in	2106395-002	
		mix 74 cm/ 29 in, 130 cm/51 in	2106395-003	

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ Monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets

Specifications available on request

Not available in all countries





ECG Leadwires Set Grabber – AHA

Description	Connection	Length	Part Number	Quantity
3-lead	Molded	74 cm/29 in	2106390-001	1
		130 cm/51 in	2106390-002	
5-lead	Non-molded	74 cm/29 in	2106391-001	1
		130 cm/51 in	2106391-002	
	Molded	74 cm/29 in	2106389-001	
		130 cm/51 in	2106389-002	
5-lead, mix	Non-molded	74 cm/29 in, 130 cm/51 in	2106391-005	1
	Molded	74 cm/29 in, 130 cm/51 in	2106389-005	
5-lead V2-6	Non-molded	130 cm/51 in	2106393-001	1
		74 cm/29 in	2106393-003	
6-lead	Non-molded	74 cm/29 in	2106397-001	1
		130 cm/51 in	2106397-002	
6-lead, mix	Non-molded	74 cm/29 in, 130 cm/51 in	2106397-003	1

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
 - B105/B125 patient monitors
 - B20/B40 and B20i/B40i patient monitors
 - Dash™ 2000/2500/3000/4000/5000 patient monitors
 - Solar™ Monitor with CARESCAPE™ Patient Data Module (PDM)
 - ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- Specifications available on request

Not available in all countries





ECG Leadwires Grabber – IEC Compliant				
Description	Connection	Length	Part Number	Quantity
3-lead	Molded	74 cm/29 in	2106390-003	1
		130 cm/51 in	2106390-004	
5-lead	Molded	74 cm/29 in	2106389-003	1
		130 cm/51 in	2106389-004	
		mix 74 cm/29 in, 130 cm/51 in	2106389-006	
	Non-molded	74 cm/29 in	2106391-003	
		130 cm/51 in	2106391-004	
		mix 74 cm/29 in, 130 cm/51 in	2106391-006	
5-lead C2-6	Non-molded	130 cm/51 in	2106393-004	1
		74 cm/29 in	2106393-002	
6-lead	Non-molded	74 cm/29 in	2106398-001	1
		130 cm/51 in	2106398-002	
		mix 74 cm/ 29 in, 130 cm/51 in	2106398-003	

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
 - B105/B125 patient monitors
 - B20/B40 and B20i/B40i patient monitors
 - Dash™ 2000/2500/3000/4000/5000 patient monitors
 - Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
 - ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- Specifications available on request

Not available in all countries





ECG Leadwires Snap – IEC Compliant				
Description	Connection	Length	Part Number	Quantity
3-lead	Molded	74 cm/29 in	2106385-003	1
		130 cm/51 in	2106385-004	
5-lead	Molded	74 cm/29 in	2106383-003	1
		130 cm/51 in	2106383-004	
		mix 74 cm/29 in, 130 cm/51 in	2106383-006	
	Non-molded	74 cm/29 in	2106381-003	
		130 cm/51 in	2106381-004	
5-lead C2-6	Non-molded	130 cm/51 in	2106394-004	1
		74 cm/29 in	2106394-002	
6-lead	Non-molded	130 cm/51 in	2106396-002	1
		74 cm/29 in	2106396-001	

Product compatibility


- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets

Specifications available on request


Not available in all countries



ECG ELECTRODES



2106891-001



2106892-001

MRI, Cath Lab, EP, Radiology –
ECG Pre-wired Radiotranslucent Electrodes with DIN Connector

Population	Conductive Medium	Cloth	Dimensions	Part Number	Quantity
PICU, NICU	Adhesive Hydrogel	X	22.225 mm/0.875 in	2106891-001* AHA	300 (100 bags of 3)
				2106892-001* IEC	

Product compatibility

GE Healthcare ECG electrodes are tested and certified compatible with GE Healthcare ECG monitoring, telemetry and diagnostic resting equipment.

Specifications available on request

** = shelf life 18 months

* = shelf life 24 months


Not available in all countries



2106924-001



2106917-001



2106921-001



2106918-001



General Purpose – ECG Snap Electrodes


Population	Conductive Medium	Foam	Cloth	Tape	Dimensions	Part Number	Quantity
Adult	Solid Hydrogel	X			38.1 mm/1.5 in	2106924-001*	600 (20 bags of 30)
			X		47.625 mm/1.875 in	2106890-001*	
		X				2106915-001*	
				X		2104777-001*	
Adult	7% KCl Wet gel	X			38.1 mm/ 1.5 in	2106921-001**	600 (10 bags of 60)
NICU, PICU	Solid Hydrogel		X		22.225 mm x 33.34 mm/ 0.875 in x 1.3125 in	2106917-001*	
		X			25.4 mm/1 in	2106918-001*	



2106889-001



2106919-001



MRI, Cath Lab, EP, Radiology – ECG Snap Electrodes Radiotranslucent

Population	Conductive Medium	Foam	Cloth	Dimensions	Part Number	Quantity
Adult	Solid Hydrogel	X		47.625 mm/1.875 in	2106889-001*	600 (10 bags of 60)
PICU NICU	Adhesive Hydrogel		X	22.225 mm/0.875 in	2106919-001*	


** = shelf life 18 months

* = shelf life 24 months

Not available in all countries



Temperature



2104178-001



2104178-002



2104179-001

Temperature Cables				
Description	Connection	Length	Part Number	Quantity
400 Series for Reusable Probes	Single	0.5 m/1.7 ft	2104178-001	1
	Dual	0.5 m/1.7 ft	2104178-002	1
400 Series for disposable Probes	Single	3.6 m/12 ft	2104179-001	1

Not available in all countries



2107176-032



2107176-033

Temperature Probes Reusable					
Description	Connection	Matching Code on the Label	Length	Part Number	Quantity
Esophageal/rectal	Adult	1010132011GE	280 cm/9 ft	2107176-031	1
	Pediatric	1010132010GE	280 cm/9 ft	2107176-032	1
Skin	Adult	1010132023GE	280 cm/9 ft	2107176-033	1

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-PSMP and E-PT modules
- B105/125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 3000/4000/5000 patient monitors
- Solar™ 8000M with TRAM module

Specifications available on request



Not available in all countries



Invasive Blood Pressure



Invasive Blood Pressure Cables				
Description	Connection	Length	Part Number	Quantity
Argon™	Single	3.6 m/12 ft	2104166-001	1
		1.2 m/4 ft	2104166-002	
	Dual	3.6 m/12 ft	2104166-003	1
		1.2 m/4 ft	2104166-004	
Edwards™ LS TruWave™	Single	3.6 m/12 ft	2104162-001	1
		1.2 m/4 ft	2104162-002	
	Dual	3.6 m/12 ft	2104162-003	1
		1.2 m/4 ft	2104162-004	

Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-PSMP, E-PP, E-PT, E-COP(SV) and E-PiCCO™ modules
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors

Specifications available on request

Not available in all countries





ICU Medical Transpac-IV® - Dual



Utah/Biotrans™ - Single



Utah/Biotrans™ - Dual



Dual Adapter Cable

Invasive Blood Pressure Cables

Description	Connection	Length	Part Number	Quantity
ICU Medical Transpac-IV®	Single	3.6 m/12 ft	2104158-001	1
		1.2 m/4 ft	2104158-002	
	Dual	3.6 m/12 ft	2104158-003	1
		1.2 m/4 ft	2104158-004	
Utah/Biotrans™	Single	3.6 m/12 ft	2104107-001	1
		1.2 m/4 ft	2104107-002	
	Dual	3.6 m/12 ft	2104107-003	1
		1.2 m/4 ft	2104107-004	
Dual Adapter Cable	–	30 cm/1 ft	2103966-001	1

Product compatibility


- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-PSMP, E-PP, E-PT, E-COP(SV) and E-PiCCO™ modules
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors

Specifications available on request

Not available in all countries



Cardiac Output



CARDIAC OUTPUT



Cardiac Output Cables				
Description	Length	For use with	Part Number	Quantity
Cardiac Output Cable	3.6 m/12 ft	CARESCAPE™ Monitors B450, B650 and B850 with the CARESCAPE™ Patient Data Module (PDM)	2104180-001	1
	2.4 m/8 ft		2104180-002	1

Specifications available on request

Not available in all countries



9446-909



9446-910

CareFusion Cardiac Output Probes					
Description	Length	For use with	Manufacturer Reference	Part Number	Quantity
Injectate Bath Probe	1.8 m/6 ft	CARESCAPE™ PDM	9446-909	9446-909	1
In-Line Injectate Probe	.9 m/3 ft	Edwards™ CO Set, PDM	9446-910	9446-910	1
	.9 m/3 ft	BD CritiKit, PDM	9446-911	9446-911	1
	2.4 m/8 ft	Thermoset, PDM	9446-913	9446-913	1

For sale with equipment purchase ONLY.

Not available in all countries

CARDIAC OUTPUT TEMPERATURE PROBES

Argon™ Accessories				
Description	For use with	Manufacturer Reference	Part Number	Quantity
Spectramed CO Set Temperature Probe	E-COP, E-COPSv	684056	16593-HEL	1

These accessories are sold ONLY with the equipment

Edwards™ Accessories				
Description	For use with	Manufacturer Reference	Part Number	Quantity
Injectate Bath Temperature Probe	E-COP, E-COPSv	9850A	16591	1
Flow-through Injectate Temperature Probe	Edwards™ catheters, E-COP, E-COPSv	93522	16592VS	1

These accessories are sold ONLY with the equipment

OTHER ACCESSORIES

Other Accessories				
Description	For use with	Manufacturer Reference	Part Number	Quantity
CO Set for Room Temperature Injectate, disposable	E-COP, E-COPSv	93610	16574VS	1
Catheter Connecting Cable	E-COP, E-COPSv	COM2-CC	16590-HEL	1
Thermodilution Catheter, disposable	E-COP, E-COPSv	131HF7	M1238258	5

These accessories are sold ONLY with the equipment

Not available in all countries




PiCCO™ Cardiac Output Cables

Description	Length	For use with	Part Number	Quantity
Continuous Cardiac Output Cable	4.1m	E-PiCCO module	M1154585	1
Injectate Sensor Cable	0.9m		M1182448	1

These accessories are sold ONLY with the equipment



M1182448



PiCCO™ Kits

Patient	Description	Includes	PULSIOCATH Arterial Thermodilation Catheters Reference	Part Number	Quantity
Adult	Femoral artery	Five Catheters, one pressure transducer with flush device and one PV4046 PiCCO™ -Technology injectate temperature sensor housing	PV2015L20 Standard	M1219758	1
	Brachial artery		PV2014L22 Brachial	M1219759	1
Children/Infant	Femoral artery		PV2013L07 Femoral	M1219760	1

These accessories are sold ONLY with the equipment

Not available in all countries



Pulse Oximetry



TruSignal™



TS-F2-GE



TS-SA4-GE



TS-E2-GE

TruSignal™ Integrated Reusable Sensors

Description	Connector	Patient	Patient Weight	Length	Part Number	Quantity
Finger	GE	Adult, Pediatric	≥ 20 kg	2 m/6.6 ft	TS-F2-GE	1
				4 m/13.1 ft	TS-F4-GE	1
	Ohmeda			1 m/3.3 ft	TS-F1-H	1
	Datex			4 m/13.1 ft	TS-F4-N	1
FingerTip	GE	Adult	> 30 kg	4 m/13.1 ft	TS-SA4-GE	1
PediTip		Pediatric	15-30 kg	3 m/9.8 ft	TS-SP3-GE	1
Ear	GE	Adult, Pediatric	≥ 10 kg	2 m/6.6 ft	TS-E2-GE	1
				4 m/13.1 ft	TS-E4-GE	1
	Ohmeda				TS-E4-H	1
	Datex				TS-E4-N	1

Not available in all countries





TS-F-D



TS-SA-D/TS-SP-D




TS-E-D



TS-W-D



TS-SE-3



TruSignal™ Reusable Sensors					
Description	Patient	Patient Weight	Length	Part Number	Quantity
Finger	Adult, Pediatric	≥ 20 kg	1 m/3.3 ft	TS-F-D	1
PediTip	Pediatric	15-30 kg		TS-SP-D	1
Fingertip	Adult	> 30 kg		TS-SA-D	
Ear	Adult, Pediatric	≥ 10 kg		TS-E-D	1
Wrap	Adult, Pediatric, Infant	≥ 3 kg		TS-W-D	1
Sensitive Skin	Adult, Pediatric, Infant, Neonate	No limit		TS-SE-3	3

Requires Compatible Interconnect Cable



TS-G3

TruSignal™ Interconnect Cables			
Connector	Length	Part Number	Quantity
GE	3 m/9.8 ft	TS-G3	1
Ohmeda		TS-H3	
Datex		TS-N3	

Not available in all countries





TruSignal™ Disposable Sensors

Description	Patient	Patient Weight	Length	Part Number	Quantity
Adult/Pediatric	Pediatric, infant	≥ 20 kg	0.3 m/ 1 ft	TS-AP-10	10
				TS-AP-25	25
AllFit	Adult, pediatric, infant, neonate	No limit	0.5 m/ 1.6 ft	TS-AF-10	10
				TS-AF-25	25
Adhesive Wrap	Adult	> 20 kg	0.3 m/ 1 ft	TS-AAW-10	10
				TS-AAW-25	25
	Pediatric, infant	3-20 kg		TS-PAW-10	10
				TS-PAW-25	25

Requires Compatible Interconnect Cable

TruSignal™ Sensor Accessories

Description	Use With	Size	Part Number	Quantity
Replacement wide adhesive tape	TruSignal™ Sensitive Skin Sensor, TruSignal™ Wrap Sensor	–	OXY-RTW	100
Replacement adhesive tape (blue)	TruSignal™ AllFit Sensor	–	OXY-RT	100
Replacement adhesive tape (teddy bears)			OXY-RTB	
Replacement foam wrap	TruSignal™ Sensitive Skin Sensor	Small	OXY-RWS	24
		Medium	OXY-RWM	
	TruSignal™ Wrap Sensor	Large	OXY-RWL	
Infant foam sandal	TruSignal™ Sensitive Skin Sensor	–	OXY-SND	3
Replacement Headband	TruSignal™ Ear Sensor	–	OXY-HB	3

Not available in all countries



For CARESCAPE™ B450 / B650 / B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Masimo™ SpO2 technology.



LNCS Digit Sensor



DCI-P Digit sensor



LNCS Tip-Clip Ear Sensor

Masimo™ Reusable LNCS and RD SET Sensors

Sensors	Masimo™ Reference	Type	Patient	Patient Weight	Part Number	Quantity
LNCS Digit Sensor	1863	DC-I	Adult	> 30 kg	207258-001	1
	1864	DCI-P	Pediatric	10-50 kg	207259-001	
LNCS Tip-Clip Ear Sensor	1895	TC-I	Adult	> 30 kg	207261-001	1
RD SET SpO2 Sensor	4050	DCI	Adult	> 30 kg	MASIMO-RD-4050	1
	4051	DCI-P	Pediatric	10-50 kg	MASIMO-RD-4051	

Requires Compatible Interconnect Cable

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.



207263-002

Masimo™ LNCS & RD SET Interconnect Cables

Description		Size	Masimo™ Reference	Part Number	Quantity
LNC-10 Interconnect Cable	For GE DINAMAP	3 m / 10 ft	1814	207263-001	1
	For GE		2016	207263-002	1
RED Interconnect Cable	LNC-04	1.5 m 5 ft	2055	2059005-001	1
	LNC-10	3 m / 10 ft	2056	2059006-001	1
	LNC-14	4.2 m / 14 ft	2057	2059007-001	1
RD SET Interconnect Cable	GE-05 (For GE)	1.5 m / 5 ft	4084	MASIMO-RD-4084	1
	GE-12 (For GE)	3.6 m / 12 ft	4085	MASIMO-RD-4085	1
RD SET Interconnect Cable	MD14-05, 14-pin	1.5 m / 5ft	4080	MASIMO-RD-4080	1
	MD14-12, 14-pin	3.6 m / 12 ft	4081	MASIMO-RD-4081	1
	MD20-05, 20-pin	1.5 m / 5 ft	4103	MASIMO-RD-4103	1
	MD20-12, 20-pin	3.6 m / 12 ft	4104	MASIMO-RD-4104	1

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

Not available in all countries





2051070-001

Masimo™ LNCs & RD SET Disposable Sensors

Description		Patient Weight	Masimo™ Reference	Part Number	Quantity
LNCs SpO2	Adult	>30 kg	1859	2027253-001	20
	Pediatric	10-50 kg	1860	2027254-001	20
	Infant	3-20 kg	2328	2051068-001	20
	Neonatal PT	< 1 kg	2330	2051070-001	20
	Adult / Neonatal	> 40 kg / < 3 kg	2329	2051069-001	20
RD SET® SpO2	Adult	>30 kg	4000	MASIMO-RD-4000	20
	Pediatric	10-50 kg	4001	MASIMO-RD-4001	20
	Infant	3-20 kg	4002	MASIMO-RD-4002	20
	Adult / Neonatal	> 40 kg / < 3 kg	4003	MASIMO-RD-4003	20
	Premature Baby	< 1 kg	4004	MASIMO-RD-4004	20

Requires Compatible Interconnect Cable

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

Masimo™ LNCs & RD SET Sensor Accessories

Description		Masimo™ Reference	Part Number	Quantity
LNCs Replacement Tape	Infant	2307	2051071-001	100
	Neonatal	2308	2051073-001	100
	Neonatal PT	2309	2051074-001	100
RD Set Replacement Tape	Infant	4046	MASIMO-RD-4046	102
	Neonatal	4047	MASIMO-RD-4047	102
RD Set Replacement Wrap	Neonatal PT	4048	MASIMO-RD-4048	10

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

Not available in all countries



For CARESCAPE™ B450 / B650 / B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Nellcor™ SpO2 technology.



407705-006



70124035



70124033



414248-001

Nellcor™ Reusable Sensors

Sensors	Patient	Patient Weight	Nellcor™ Reference	Part Number	Quantity
DuraSensor Finger Probe	Adult	< 30 kg	DS100A-1	407705-006	1
Multisite SpO2 Sensor	Adult	> 1 kg	D-YS	70124033	1
Two-piece SpO2 sensor	Neonatal / Adult	< 3 kg / > 40 kg	OXI-A/N	70124035	1 (delivered with 24 adhesive strips)
	Pediatric, infant	3-40 kg	OXI-P/I	414248-001	1 (delivered with 24 adhesive strips)

Requires Compatible Interconnect Cable

These accessories are sold ONLY with the equipment

Not available in all countries





Nellcor™ Interconnect Cables				
Description	Connector	Length	Part Number	Quantity
OxiMax™ SpO2 Interconnect Cable	GE	1.2 m / 4 ft	2021406-002	1
		3 m / 10 ft	2021406-001	1
	GE for Corometrics™ only	2.9 m / 10 ft	2025350-001	1
OxiSmart SpO2 Interconnect Cable	GE	2.9 m / 10 ft	2006644-001	1
Interconnect Cables*	DOC-4	1.2 m / 4 ft	2059000-001	1
	DOC-8	2.4 m / 8 ft	2059001-001	1
	DOC-10	3 m / 10ft	2008773-001	1

* These accessories are sold ONLY with the equipment

Not available in all countries

For CARESCAPE™ B450 / B650 / B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Nellcor™ SpO2 technology.



70124027



70124032

Nellcor™ Disposable Sensors

Sensors	Nellcor™ Reference	Patient	Patient Weight	Part Number	Quantity
Adhesive SpO2 sensor	MAX-A	Adult	> 30 kg	70124027	24
Long Finger Adhesive sensor	MAXAL			2028117-001	24
Adhesive SpO2 sensor	MAXP	Pediatric	10-50 kg	70124022	24
	MAXI	Infant	3-20 kg	70124026	24
	MAX-N	Adult	> 40 kg	70124032	24
		Neonate	< 3 kg		

Requires Compatible Interconnect Cable
These accessories are sold ONLY with the equipment



Nellcor™ Fixation Clip

Description	Nellcor™ Reference	For Use With	Part Number	Quantity
SpO2 Ear-Clip	D-YSE	D-YS	70124034	1

These accessories are sold ONLY with the equipment



Nellcor™ Sensor Adhesive Tape

Product	Patient	Nellcor™ Reference	Use With	Part Number	Quantity
Sensor Wrap Replacement Tape	Adult/Neonatal	ADH-A/N	OXI-A/N, D-YS	2016130-001	100
	Pediatric/Infant	ADH-P/I	OXI-P/I, D-YS	2016131-001	100

These accessories are sold ONLY with the equipment

Not available in all countries



Adequacy of Anesthesia



ENTROPY



M1038681

M1174413

Sensors			
Product	Use	Part Number	Quantity
Entropy™ EasyFit Sensor	Disposable	M1174413	25
Entropy™ Sensor		M1038681	



Cable			
Product	Length	Part Number	Quantity
Entropy™ Cable	3.5 m	M1050784	1

Product compatibility

E-ENTROPY module

The GE Healthcare Entropy™ module, E-ENTROPY, and accessories are indicated for adult and pediatric patients older than 2 years within a hospital for monitoring the state of the brain by data acquisition of electroencephalograph (EEG) and frontal electromyograph (FEMG) signals during general anesthesia.

Not available in all countries



Electromyography (EMG)
Requires 5 electrodes



Kinemyography (KMG)
Requires 2 electrodes



888418/897439



888416

Sensors				
Product	Use	Patient	Part Number	Quantity
NMT MechanoSensor	Reusable	Adult	888418	1
		Pediatric	897439	
NMT ElectroSensor		Adult, Pediatric	888416	1
NMT Regional Block Adapter			888417	1

Accessories				
Product	Use	Patient	Part Number	Quantity
NMT Electrodes	Disposable	Adult, Pediatric	57268-HEL	30

Cable			
Product	Length	Part Number	Quantity
NMT Sensor Cable	3 m	888414	1

Product compatibility

E-NMT module

Not available in all countries



Neurology



Reusable Preconfigured EEG Leadsets			
Product	Length	Part Number	Quantity
Basic montage	0.6 m	898050	1
General montage		898051	1
AEP montage		898052	1

Not available in all countries





2083319-001



2083320-001

Disposable EEG Electrodes

Product	Part Number	Quantity
EEG Clip electrodes Ag/AgCl	M1240794	40 pouches of 25 (1,000 total)
	572686	500
EEG Prewired cup electrodes	2083319-001	46 pouches of 10 (460 total)
EEG Prewired wet gel electrodes	2083320-001	42 pouches of 10 (420 total)

All electrodes within a set-up should be similar

Product compatibility

E-EEG(X) module and N-EEG(X) headbox

AEP Accessories

Product	Part Number	Quantity
AEP Earphone set including earphone, sound divider, cable package and disposable ear tips	8001000	1 earphone 20 disposable ear tips
Disposable ear tips for use with earphone set	8001233	20

To use with 898052

Product compatibility

E-EEG(X) module and N-EEG(X) headbox

Not available in all countries



Patient Spirometry



896952



8001948



733910-HEL



73393

Sensors				
Description	Type	Approved with Modules	Part Number	Quantity
D-lite™ +	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE¹, E-sCAiOVX	896952	50
Pedi-lite+			8001948	
D-lite™	Re-usable		733910-HEL	1
Pedi-lite			73393	

Specifications available on request

(1) Not 510(k) cleared

Not available in all countries





Spirometry Tubes

Description	Length	Type	Approved with Modules	Part Number	Quantity
Spirometry Tubes	2m/7ft	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE ¹ , E-sCAiOVX, E-COV, E-COVX, E-CAiOV, E-CAiOVX	890031	5
	3m/10ft			884101	
	6m/20ft			891191-HEL	



Kits

(include gas sample line, sensor and spirometry tube)

Description	Length	Type	Approved with Modules	Part Number	Quantity
D-lite™ ++ kit	2m/7ft	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE ¹ , E-sCAiOVX	2104297-001	20
	3m/10ft			2104297-002	
D-lite™ + kit	2m/7ft	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE ¹ , E-sCAiOVX, E-COV, E-COVX, E-CAiOV, E-CAiOVX	8004381	20
	3m/10ft			8004382	
Pedi-lite+ kit	2m/7ft	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE ¹ , E-sCAiOVX, E-COV, E-COVX, E-CAiOV, E-CAiOVX	8002718	20
	3m/10ft			M1032634	



Specifications available on request

(1) Not 510(k) cleared

Not available in all countries

Gas Monitoring



ICU AIRWAY ADAPTER SETS



Airway Gas Sampling Sets

Description	Type	Approved with Modules	Length	Part Number	Quantity
Adult/Pediatric	Single Patient Disposable	E-miniC, E-sCO, E-sCOV, E-sCOVX	3 m/10 ft	8090121531	20 pcs/box
			2 m/7 ft	8090121521	
Infant		E-sCO, E-sCOV, E-sCOVX	2 m/7 ft	8090121522	20 pcs/box

Not available in all countries



GAS SAMPLING LINES AND ADAPTERS



Gas Sampling Lines					
Description	Type	Approved with Modules	Length	Part Number	Quantity
PVC/PE gas sampling line	Single Patient Disposable	E-miniC, E-sCO, E-sCOV, E-sCOVX, E-sCAiO, E-sCAiOE ¹ , E-sCAiOV, E-sCAiOVX, E-sCAiOVE ¹ , N-CAiO	2 m/7 ft	2097307-001	10
			3 m/7 ft	2097307-002	
			6 m/20 ft	2097307-003	5



Nasal Gas Sampling Lines – Vyaire					
Description	Type	Approved with Modules	Length	Part Number	Quantity
Pediatric	Single Patient Disposable	E-miniC, E-sCO, E-sCOV, E-sCOVX, E-sCAiO, N-CAiO E-sCAiOE ¹ , E-sCAiOV, E-sCAiOVX, E-sCAiOVE ¹	2 m/7 ft	878033	10
Adult				878034	
Large Adult				878035	

These accessories are sold ONLY with the equipment



73385-HEL

Gas Monitoring Airway Adapters – Vyaire				
Description	Type	Approved with Modules	Part Number	Quantity
Straight T-adapter	Single Patient Disposable	E-miniC, E-sCO, E-sCOV, E-sCOVX, E-sCAiO, E-sCAiOE ¹ , E-sCAiOV, E-sCAiOVX, E-sCAiOVE ¹	73385-HEL	10
Elbow adapter			73386-HEL	

These accessories are sold ONLY with the equipment



3090121002

Gas Monitoring Airway Adapter				
Description	Type	Approved with Modules	Part Number	Quantity
Airway Adapter with Luer connector	Single Patient Disposable	E-miniC, E-sCO, E-sCOV, E-sCOVX, E-sCAiO, E-sCAiOE ¹ , E-sCAiOV, E-sCAiOVX, E-sCAiOVE ¹ , N-CAiO	3090121002	25

(1) Not 510(k) cleared

Not available in all countries



NEONATAL LOW DEAD SPACE ADAPTERS



8001540

Neonatal Low Dead Space Adapters – Vyaire						
Description	ET Tube ID	Dead Space	Type	Approved with Modules	Part Number	Quantity
Neonatal Low Dead Space Adapters, Single Patient Disposable, 2.5mm	3 mm	0.15 ml	Single Patient Disposable	E-sCO, E-sCOV, E-sCOVX, E-sCAiO, E-sCAiOE ¹ , E-sCAiOV, E-sCAiOVX, E-sCAiOVE ¹ , N-CAiO	8001540	10
	3 mm	0.15 ml		E-miniC, E-sCO, E-sCOV, E-sCOVX, E-sCAiO, E-sCAiOE ¹ , E-sCAiOV, E-sCAiOVX, E-sCAiOVE ¹ , N-CAiO	8001541	10
	3.5 mm	0.2 ml			8001542	
	4 mm	1 ml			8001543	
	4.5 mm				8000062	
	5 mm				8000063	

These accessories are sold ONLY with the equipment

WATER TRAPS



881319-HEL



M1200227



876446-HEL



M1182629



8002174

Water Traps – ICU

Description	Type	Approved with Modules	Part Number	Quantity
D-fend™ +	Single Patient Disposable	E-CO, E-COV, E-COVX	881319-HEL	10
D-fend™ Pro+		E-sCO, E-sCOV, E-sCOVX	M1200227	
Mini D-fend™		E-miniC	8002174	

Specifications available on request

Water Traps – Anesthesia

Description	Type	Approved with Modules	Part Number	Quantity
D-fend™	Multi-Patient Disposables	E-CAiO, E-CAiOV, E-CAiOVX	876446-HEL	10
D-fend™ Pro		E-sCAiO, E-sCAiOE ¹ , E-sCAiOV, E-sCAiOVX, E-sCAiOVE ¹	M1182629	
Mini D-fend™		E-miniC	8002174	

Specifications available on request

(1) Not 510(k) cleared

Not available in all countries



GAS SCAVENGING AND REGULATORS



M1031274



8004462

Gas Scavenging – Vyaire

Description	Type	Approved with Modules	Length	Part Number	Quantity
Exhaust Line for gas return or scavenging	Single Patient Disposable	E-miniC, E-sCAiO, E-sCAiOE ¹ , E-sCAiOV, E-sCAiOVX, E-sCAiOVE ¹	2 m/7 ft	M1031274	1
Exhaust Line with Colder Fitting		E-sCAiO, E-sCAiOE ¹ , E-sCAiOV, E-sCAiOVX, E-sCAiOVE ¹	1 m/41 in	8004462	1
			18 cm/7 in	8004463	

These accessories are sold ONLY with the equipment

Calibration Gases and Regulators

Description	CO2	N2O	O2	Desflurane	N2	Type	Approved with Modules	Part Number	Quantity
Calibration Gas, (with Regulator 755534-HEL)	5%	33%	55%	2%	5%	Disposable	E-sCAiO, E-sCAiOV, E-sCAiOVX, E-sCAiOE ¹ , E-sCAiOVE ¹ , N-CAiO, E-CAiO, E-CAiOV, E-CAiOVX	755583-HEL**	1
Calibration Gas, (with Regulator 755534-HEL)		40%	55%	–	–		E-sCO, E-sCOV, E-sCOVX, E-CO, E-COV, E-COVX	755581-HEL	4
Quick Calibration Gas, (with Regulator 755534-HEL)		75.05%	19.95%	–	–		E-miniC	755580	4
Calibration Gas		36%	54.5%	2%	–		E-sCAiO, E-sCAiOV, E-sCAiOVX, N-CAiO, E-CAiOVX, E-CAiO, E-CAiOV	755571-HEL*	1
Quick Calibration Gas, (With Regulator M1006864)		–	95%	–	–		E-sCO, E-sCOV, E-sCOVX, E-CO, E-COV, E-COVX	755587*	4

*Available in US only

**Not available in US

(1) Not 510(k) cleared

Not available in all countries



ANESTHESIA & RESPIRATORY CARE



ANESTHESIA & RESPIRATORY CARE

Absorbents

- Absorbers
- Canisters
- Dust Filter

Vaporizer Bottle Adapters

Filters & Breathing Masks

Reservoir Bags & Gas Scavenging

- Anesthesia Breathing Bags
- Scavenging Bags, Tubes, & Connectors

Patient Breathing Circuits

- Circuits
- Circuit Components
- Circuit Kits
- Circuit Startup Kits

Flow Sensors - Respiratory

- Cables
- Heater
- Respiratory Flow Sensors and Expiratory Valve Assemblies

Flow Sensors - Anesthesia

Nebulizer Accessories

Anesthesia O2 Cells

Respiratory Safety Guard

Pressure Catheters

Anesthesia & Respiratory care



Disposable CO2 Absorbers

Description	Compatibility	Part Number	Quantity
AMSORB® Plus G-Can Absorber, 1.0L	All Aisys™, Avance™, Aespire™, 9100c and ADU variants	2105489-003	8
AMSORB® Plus Prefilled Cartridge, 1.2L	All Aestiva™ variants	2105489-004	12
AMSORB® Plus CareCan Absorber, 1.4L	Carestation™ 620 & 650	2105489-006	8

Reusable CO2 Absorbent Components

Description	Compatibility	Part Number	Quantity
AMSORB® Plus Jerican Loosefill Absorbent, 5.0L	GE Healthcare Reusable Canisters	2105489-007	2
ABS Reusable CO2 Canister	All Aisys™, Avance™, Aespire™ & 9100c variants	1407-7004-000	1
CBS Reusable CO2 Canister	Carestation™ 620 & 650	2071165-001-S	1
Foam Dust Filter	ABS Reusable Canister	2105489-001	40

Easy-Fil Vaporizer Bottle Adapters

Vaporizer Compatibility	Color Indicator	Part Number	Quantity
Isoflurane Easy-Fil	Purple	2106572-001	1
Halothane Easy-Fil	Red	2106572-002	
Sevoflurane Easy-Fil	Yellow	2106572-003	





HMEFs and Filters – Disposable

Description	Size	Sampling Port	Part Number	Quantity
HMEF	Neonatal	Included	2106570-006	50
	Pediatric		2106570-009	
	Adult		2106570-010	
Bacterial and Viral Filter	Adult	Included	2106570-007	50
HEPA Filter	Adult		2106570-008	



Anesthesia Breathing Masks – Disposable

Mask Connection	Mask Size	Part Number	Quantity
15mm	0 – Neonatal	2106570-012	50
	1 – Small Child	2106570-013	
22mm	2 – Child	2106570-014	50
	3 – Large Child	2106570-015	
	4 – Adult	2106570-016	50
	5 – Large Adult	2106570-017	
	6 – Extra Large Adult	2106570-018	



Anesthesia Breathing Masks – Reusable

Mask Connection	Mask Size	Part Number	Quantity
15mm	0 – Neonatal	2106570-011	1
		2106570-029	10
	1 – Small Child	2106570-019	1
		2106570-030	10
22mm	2 – Child	2106570-020	1
		2106570-031	10
	3 – Large Child	2106570-021	1
		2106570-032	10
	4 – Adult	2106570-022	1
		2106570-033	10
	5 – Large Adult	2106570-023	1
		2106570-034	10

Not available in all countries



RESERVOIR BAGS AND GAS SCAVENGING



Anesthesia Manual Breathing Bags – Reusable				
Bag Connection Size	Material	Bag Size	Part Number	Quantity
22mm	Clear Silicone	0.5 L	2106570-001	1
			2106570-025	5
		1.0 L	2106570-002	1
			2106570-026	5
		2.0 L	2106570-003	1
			2106570-027	5
		3.0 L	2106570-004	1
			2106570-028	5



Passive Scavenging Tubes – Disposable				
Connection	Length	Type	Part Number	Quantity
30M-30M	5 m/ 16.4 ft	Tube	2105488-006	1
			2105488-022	10



Gas Scavenging Bag – Reusable				
Size	Compatibility	Connection	Part Number	Quantity
3.0L	Aisys, Avance, Aespire™ (ABS Breathing Systems)	30M	2106570-024	1

Straight Connectors – Scavenging			
Connection	Reprocessable	Part Number	Quantity
30M-19M	Disposable	M1003947	1
		1500-3376-000	5
30M-30M		M1003134	1

Not available in all countries



PATIENT BREATHING CIRCUIT KITS – START UP

All startup kits include patient breathing circuit, reservoir bag tube, reservoir bag, gas sampling lines and gas sampling exhaust lines. These are primarily intended for initial anesthesia machine setup.



2096534-009
Reusable, Straight Adapter



2096536-002
Disposable, non-Expandable, ABS Adapter



2105488-019
Disposable, Single Limb, ABS Adapter



2096533-001
Disposable, ExPanda™able, Straight Adapter

Patient Breathing Circuit Kits – Start Up – Adult

Type	Gas Sampling Line	Gas Sampling Exhaust Line	Limb Circuit	Breathing Circuit	Breathing Bag		Part Number	Quantity
				Length	Bag Connector	Size		
Disposable	3m PVC/PE	0.18m & 1m	Dual	Expandable up to 3 m/9.8 ft	Straight	3L	2096533-001	1
						3L	2096533-002*	1
						2L	2096533-003	1
			Single Split Lumen	1.8 m/5.9 ft	Angled (ABS)	2L	2096536-002	2
					Straight	2L	2096536-003	2
					Straight	3L	2096536-004*	2
					Angled (ABS)	2L	2105488-019	1
Reusable			Dual	1.5 m/4.9 ft	Angled (ABS)	2L	2105488-024	20
					Straight	2L	2096534-009	1
					Angled (ABS)	2L	2096536-001	1

*Australia and New Zealand ONLY

Not available in all countries



PATIENT BREATHING CIRCUITS

Our wide range of patient breathing circuits and components offer robust and high quality options for both disposable and reusable applications.

Patient Breathing Circuit Kits									
Type	Size	3m PVC/PE Gas Sampling Line	Circuit Type	Breathing Circuit	Breathing Bag		Includes	Part Number	Qty
				Length	Bag Connector	Size			
Disposable	Adult	Yes	Dual Limb	Expandable up to 3 m/9.8 ft	Not included	3L	-	2096535-006	20
				1.8 m/5.9 ft	Angled (ABS)		-	2105488-013	20
					Straight		-	2105488-014	20
		No	Single Split Lumen	1.8 m/5.9 ft	Straight	3L	-	2105488-010	20
			Dual Limb	Expandable up to 1.8 m/5.9 ft	Not included		-	2105488-011	20
				1.5 m/4.9 ft	Not included		2x Water Traps	2096535-001	20
							1x Water Trap	2105488-017	20
			Not Included		Angled (ABS)	2L	Bag Tube & Bag Only	2105488-012	20
	Neonatal	Yes	Dual Limb	1.5 m/4.9 ft	Not included	1L	HME & Filter	2105488-001	20
		No						2105488-021	1
Reusable	Adult	No	Dual Limb	1.5 m/4.9 ft	Not included		2x Water Traps	2096534-005	1
	Pediatric							2096534-006	1



Not available in all countries



PATIENT BREATHING CIRCUITS



2096534-007
Reusable, Pediatric



2105488-020
Disposable, Pediatric



2105488-015
Disposable, Adult



2105488-009
Disposable, Adult, Inteliflo® Single Split Lumen

Patient Breathing Circuits

Type	Size	Circuit Type	Breathing Circuit Length	Patient Elbow	Gas Sampling Port Location	Part Number	Quantity
Reusable	Adult	Dual Limb	1.5 m/ 4.9 ft	Double Swivel, Removable	Elbow	2096534-007	1
	Pediatric					2096534-008	1
Disposable	Adult	Single Split Lumen	1.8 m/ 5.9 ft	Fixed, Removable	Elbow	2105488-009	20
		Dual Limb	1.5 m/ 4.9 ft	Double Swivel, Removable	Patient Wye	2096535-002	20
				Fixed, Removable	Elbow	2105488-015	20
	Pediatric			Double Swivel, Removable	Patient Wye	2105488-020	20



Non-Recirculating Breathing Circuits – Jackson Rees


Type	Breathing Circuit Length	Bag Size	Gas Sampling Port Location	Part Number	Quantity
Disposable	0.8 m/2.6 ft	1 L	Elbow	2105488-008	20

Not available in all countries



PATIENT BREATHING CIRCUITS

Reusable Patient Breathing Circuit Tubes					
Size	Material	Length	Connections	Part Number	Quantity
Adult	Hytrel®	0.35 m/1.2 ft	22F-22F (90 Deg)	2096534-001	1
		0.75 m/2.4 ft		2096534-002	1
		0.9 m/3 ft		2096534-003	1
		1.5 m/4.9 ft		2096534-004	1
	Silicone	1 m/3.3 ft	22F-22F	2096534-101	1
		1.5 m/4.9 ft		2096534-102	1
Pediatric	Silicone	1.5 m/4.9 ft	22F-10F	2096534-103	1
	Hytrel®	1.1 m/3.6 ft	22F-15F	2105488-007	1
				2105488-023	10



Not available in all countries



Patient Breathing Circuit Components				
Description	Size	Reprocessable	Part Number	Quantity
Y-Piece Connector Kit, 15M-10F	Neonatal	Reusable	2096534-104	10
Patient Angled Y-Piece			2105488-002	10
Patient Straight Y-Piece - Swivel			2105488-004	10
Straight Connector - CBS Bag Arm, 22M-22M	Adult	Reusable	2096533-004	1
ABS Bag Arm Connector, 22M-22M	All		2105488-005	10
Patient Elbow - Swivel	Adult	Reusable	2105488-003	10
Patient Elbow		Disposable	2105488-018	50
Gas Return Adapter, 22M-22F	Adult	Disposable	2105488-016	10

Not available in all countries

FLOW SENSORS – RESPIRATORY



Respiratory Flow Sensors and Expiratory Valve Assemblies (EVA)				
Description	Reprocessable	Patient Type	Part Number	Quantity
Flow Sensor Only	Reusable	All	1505-3231-000	1
EVA + Flow Sensor			1505-3848-000	1
EVA Only			1505-8568-000	1
Flow Sensor Only (Neonatal)		Neonatal	1505-3272-000	1
Flow Sensor Only	Single Patient Use	All	1505-3231-001	1
EVA + Flow Sensor			1505-3848-001	1
EVA Only			1505-8568-001	1
Flow Sensor Only			1505-3231-010	10
EVA + Flow Sensor			1505-3848-010	10
EVA Only			1505-8568-010	10

Product compatibility
CARESCAPE™ R860 and Engström Ventilators

Not available in all countries




FLOW SENSORS – RESPIRATORY



Respiratory Flow Sensor and Expiratory Valve Components			
Description	Patient Type	Part Number	Quantity
Flow Sensor Cable	Neonatal	1505-5604-000	1
Exhalation Valve Heater	All	M1200693	1
Exhalation Valve Heater Cable		M1188723	1

Product compatibility
 CARESCAPE™ R860 and Engström Ventilators

FLOW SENSORS – ANESTHESIA



Anesthesia Flow Sensors			
Machine Compatability	Autoclavable	Part Number	Quantity
Aisys™, Avance™, Aespire™, Aestiva™	Yes	2087640-001-S	1
	No	2089610-001-S	
Carestation™ 600 Series	Yes	2069358-001-S	1

Not available in all countries

NEBULIZER ACCESSORIES



Aerogen® Nebulizer Accessories			
Description	Patient Type	Part Number	Quantity
Aerogen® Pro Nebulizer kit with adult and ped tee and cable	Adult and Pediatric	1505-3846-000	1
Nebulizer head with filler cap	All	AG-AP1000	
Nebulizer Cable		1505-5602-000	
Nebulizer tee with plug	Adult	AG-AP1010	5
	Pediatric	AG-AP1020	
	Neonatal	AG-AP1035	
Nebulizer replacement kit	Adult	AG-AP1100	1
Aerogen® Solo Nebulizer		AG-AS3100-S	5
		AG-AS3200-S	10
Nebulizer tee	Adult	AG-AS3010	10
	Pediatric	AG-AS3020	
	Neonatal	AG-AS3035	

Product compatibility

CARESCAPE™ R860 and Engström Ventilators

These accessories are sold ONLY with the equipment.

ANESTHESIA O2 CELLS

Anesthesia O2 Cells			
Description	Machine Compatibility	Part Number	Quantity
O2 Sensor ONLY	Aisys, Avance, Aespire™	M-10	1
O2 Sensor and cable		2082006-001	
O2 Sensor with Adapter and O-Ring	Carestation™ 600 Series	2071018-001-S	



2071018-001-S

Not available in all countries



RESPIRATORY INSPIRATORY SAFETY GUARD



Respiratory Inspiratory Safety Guard			
Machine Compatability	Patient Type	Part Number	Quantity
CARESCAPE™ R860 and Engström Ventilators	All	2066713-001	1
		2083208-001	10

PRESSURE CATHETERS




Respiratory Pressure Catheters			
Description	Machine Compatability	Part Number	Quantity
Intratracheal Pressure Catheter	CARESCAPE™ R860 and Engström Ventilators	M1045564	10

Used in conjunction with Spirodynamics Software.

Not available in all countries





DIAGNOSTIC CARDIOLOGY

MAC™, CASE™ and CardioSoft™

ECG

- Accessories
- Cables
- Electrodes
- Leadwires
- Recording Paper

SEER® Holter Monitors

ECG – CM 4000

- Leadwire Sets
- Accessories

ECG – SEER® 12

- Leadwire Sets
- Accessories


ECG – SEER® 1000

- Leadwire Sets
- Accessories

ECG – SEER® Light

- Leadwire and Cable Kits
- Accessories

MAC™, CASE™ and CardioSoft™



CABLES AND LEADWIRES – AHA



ECG Trunk Cables – AHA

Description	Length	Part Number	MAC™						CASE™	CardioSoft™	Quantity
			VU360	5500	3500	2000	800	600			
10-lead AHA	220 cm / 87 in	2104723-001			X	X	X	X			1

Integrated Cable and Leadwire Set – AHA

Description	Length	Part Number	MAC™						CASE™	CardioSoft™	Quantity
			VU360	5500	3500	2000	800	600			
ECG Value Cable 10-lead Cable / Leadwire AHA	Limb leads: 125 cm / 49 in Chest leads: 80 cm / 32 in Trunk cable: 220 cm / 87 in	2104727-001				X	X	X			1



Not available in all countries





2104724-001



2104751-001



2104752-001

ECG Leadwire Set – Banana – AHA

Description	Length	Part Number	MAC™						CASE™	CardioSoft™	Quantity
			VU360	5500	3500	2000	800	600			
Base 10, banana, AHA	Limb leads: 125 cm / 49 in Chest leads: 80 cm / 32 in	2104724-001	X	X	X	X	X	X	X	X	1
4 Add-on, banana, AHA-HEIM		2104751-001	X	X					X	X	1
3 Add-on, banana, AHA-PED (V3R, V4R, V7)		2104752-001	X	X					X	X	1



2104749-001



2104725-001

ECG Leadwire Set – Grabber – AHA

Description	Length	Part Number	MAC™						CASE™	CardioSoft™	Quantity
			VU360	5500	3500	2000	800	600			
Base 10, grabber, AHA	Limb leads: 125 cm / 49 in Chest leads: 80 cm / 32 in	2104749-001	X	X	X	X			X	X	1
4 Add-on, grabber, AHA-HEIM		2104725-001	X	X					X	X	1

Not available in all countries



CABLES AND LEADWIRES – IEC

ECG Trunk Cables – IEC

Description	Length	Part Number	MAC™						CASE™	CardioSoft™
			VU360	5500	3500	2000	800	600		
Trunk Cable 10-lead IEC	220 cm / 87 in	2104722-001			X	X	X	X		



Integrated Cable and Leadwire Set – IEC

Description	Length	Part Number	MAC™						CASE™	CardioSoft™
			VU360	5500	3500	2000	800	600		
ECG Value Cable 10-lead Cable / Leadwire IEC	Chest leads: 80 cm / 32 in	2104726-001								
	Limb leads: 125 cm / 49 in					X	X	X		
	Trunk cable: 220cm / 87 in									



ECG Leadwire Set – Banana – IEC

Description	Length	Part Number	MAC™						CASE™	CardioSoft™
			VU360	5500	3500	2000	800	600		
Leadwire Set, Base 10, banana, IEC	Limb leads: 125 cm / 49 in	2104724-002	X	X	X	X	X	X	X	X
Leadwire Set, 4 Add-on, banana, IEC-HEIM	Chest leads: 80 cm / 32 in	2104750-001	X	X					X	X



Not available in all countries



CABLES AND LEADWIRES – IEC



2104748-001



2104725-002

ECG Leadwire Set – Grabber – IEC

Description	Length	Part Number	MAC™						CASE™	CardioSoft™
			VU360	5500	3500	2000	800	600		
Base 10, IEC	Limb lead: 125 cm / 49 in	2104748-001	X	X	X	X			X	X
4 Add-on, IEC-HEIM	Chest lead: 80 cm / 32 in	2104725-002	X	X					X	X



2104729-001



2104767-001

Accessories

Description	Part Number	MAC™						CASE™	CardioSoft™	Quantity
		VU360	5500	3500	2000	800	600			
Adapter Clip, Wide Mouth Clip	2104729-001	X	X	X	X	X	X	X	X	14
ECG Leadwire Port Plug	2104767-001	X	X					X	X	4

CARDIOSOFT™



CardioSoft™ Accessories

Description	Part Number	Quantity
Supply BP Bag Belt TONOPORT V	2104824-001	1

Not available in all countries




ECG RECORDING PAPER



DCAR ECG Recording Thermal Paper													
Folding	White patient data area	Paper size	Red grid wide	Queue	MAC™						CASE™	Part Number	Quantity
					VU360	5500	3500	2000	800	600			
z-fold	No	A4 *	195 mm/ 7.6 in	Hole		X	X				X	2104769-001	300 sheets/ pack, 8 packs
			205 mm/ 8 in	Block	X			X				2104772-001	150 sheets/ pack, 10 packs
		21.5 cm x 28 cm/ 8.5 in x 11 in	200 mm/ 7.8 in	Hole		X	X				X	2104768-001	300 sheets/ pack, 8 packs
		80 mm x 90 mm/ 3.1 in x 3.5 cm	75 mm/ 2.9 in	Block						X		2104773-001	280 sheets/ pack, 10 packs
		110 mm x 140 mm/ 4.3 in x 5.5 cm	100 mm/ 3.9 in	Block					X			2104775-001	200 sheets/ pack, 10 packs
	Yes	21.5 cm x 28 cm/ 8.5 in x 11 in	155 mm/ 6.1 in	Hole	X	X	X					2104770-001	300 sheets/ pack, 8 packs
				Hole	X		X	X				2104771-001	150 sheets/ pack, 16 packs

For ECG paper for Patient monitoring and Diagnostic cardiology equipment, please go to pages 6 and 72.



GE Healthcare reusable bulb electrodes with silver-silver chloride (Ag/AgCl) plating are designed to ensure positive contact with the skin and provide maximum signal quality.



Reusable Bulb Electrode

Description	Part Number	Case Quantity
Reusable Bulb Electrode	2104783-001	6



Bulb Electrode Accessories

Description	Compatibility	Part Number	Case Quantity
Clamp Limb AHA Value	MAC™ 600, MAC™ 800, MAC™ 2000 ECG System	2104784-001	4/set
Clamp Limb IEC Value		2104785-001	

Specifications available on request

Not available in US



Stress and Holter Testing – ECG Snap Electrodes


Population	Conductive Medium	Foam	Dimensions	Part Number	Quantity
Adult	10% KCI Wet gel	x	38.1 mm/1.5 in	2104779-001**	600 (10 bags of 60)

For ECG Electrodes for Patient Monitoring, go to page 15.

Product compatibility

GE Healthcare ECG electrodes are tested and certified compatible with GE Healthcare ECG monitoring, telemetry and diagnostic resting equipment.

Specifications available on request



Diagnostic Resting – ECG Tab Electrodes

Population	Conductive Medium	Adhesive	Dimensions	Part Number	Quantity
Adult	Solid Hydrogel	x	Large: 28.575 mm x 36.525 mm/1.125 x 1.438 in	2104780-001*	4000 (4 boxes of 1000)
			Medium: 34.036 mm x 22.86 mm/ 1.34 in x .90 in	2106911-001*	
PICU NICU	Adhesive Hydrogel	x	Small: 12.7 mm x 32.54 mm/0.5 in x 1.281 in	2104778-001*	

** = shelf life 18 months

* = shelf life 24 months

Not available in all countries



SEER® Holter Monitors



CM 4000 LEADWIRE SETS



2105803-006



2105803-005

CM 4000 Leadwire Set – AHA

Description	Length	Part Number	Quantity
3-leads	80 cm/31.5 in	2105803-002	1
5-leads		2105803-004	1
7-leads		2105803-006	1

Specifications available on request

CM 4000 Leadwire Set – IEC

Description	Length	Part Number	Quantity
3-leads	80 cm/31.5 in	2105803-001	1
5-leads		2105803-003	1
7-leads		2105803-005	1

Sold ONLY in Europe

Specifications available on request

Not available in all countries



CM 4000 ACCESSORIES

CM 4000 Accessories

Description	Part Number	Case Quantity
Carry Pouch Reusable	2054933-014	1
Carry Bag Disposable	2054933-015	1
USB Cable	2054933-013	1

Specifications available on request

SEER® 12



2107228-001



2107229-001



2104808-001

SEER® 12 – AHA

Description	Length	Type	Part Number	Quantity
12-Channel Holter ECG Cable	130 cm/51.2 in	Round	2107228-001	1
	105 cm/41.3 in	Flat	2107229-001	

SEER® 12 ECG cables are compatible with SEER® 12 monitors and CardioMem™ CM 3000-12 BT monitors

SEER® 12 – IEC

Description	Length	Type	Part Number	Quantity
12-Channel Holter ECG Cable	130 cm/51 in	Round	2104808-001	1
	105 cm/41 in	Flat	2104809-001	

SEER® 12 ECG cables are compatible with SEER® 12 monitors and CardioMem™ CM 3000-12 BT monitors



2104803-001

SEER® 12 Holter Monitor Accessories

Description	Part Number	Quantity
SEER® 12 lead portable pouch, waist belt, shoulder strap, 1 battery compartment	2104803-001	1
SEER® 12 Carrying pouch, 2 battery compartment	2014606-297	1

2104803-001 is compatible with SEER® 12, CM 3000, CM 3000-12
2014606-297 must be purchased with service
Specification available on request

Not available in all countries





SEER® 1000 Leadwire Set – AHA

Description	Length	Part Number	Quantity
3 leadwire, 3 CH	75 cm/29.5 in	2104821-001	1
	105 cm/41 in	2104822-001	
5 leadwire, 2 CH	75 cm/29.5 in	2104814-001	
	105 cm/41 in	2104816-001	
7 leadwire, 3 CH	75 cm/29.5 in	2104815-001	
	105 cm/41 in	2104817-001	

Product compatibility

SEER® 1000 ECG cables are compatible with SEER® 1000 24-hour Holter monitor, SEER® 1000 48-hour Holter monitor and SEER® 1000 7-day Holter monitors

Specifications available on request

SEER® 1000 Leadwire Set – IEC

Description	Length	Part Number	Quantity
5 leadwire, 2 CH	75 cm/29.5 in	2104810-001	1
7 leadwire, 3 CH		2104811-001	
5 leadwire, 2 CH	105 cm/41 in	2104812-001	
7 leadwire, 3 CH		2104813-001	
3 leadwire, 3 CH	75 cm/29.5 in	2104818-001	
	105 cm/41 in	2104819-001	

Product compatibility

SEER® 1000 ECG cables are compatible with SEER® 1000 24-hour Holter monitor, SEER® 1000 48-hour Holter monitor and SEER® 1000 7-day Holter monitors

Specifications available on request

Not available in all countries



2067634-073

SEER® 1000 Holter Monitor Accessories		
Description	Part Number	Quantity
Pouch with waist belt and shoulder strap	2067634-073	1
Pouch Belt Clip	2067634-166	1

Specification available on request



2067634-166

Not available in all countries

SEER® Light Leadwire Set Kits – AHA			
Description	Length	Part Number	Quantity
Carry Pouch and 2 Channel Patient cable	140 cm/55 in	2042686-001	1
Carry Pouch and 3 Channel Patient cable		2042686-002	
Carry Pouch and 2 Channel Patient cable	60 cm/24 in	2042686-005	1
Carry Pouch and 3 Channel Patient cable		2042686-006	

Product compatibility

SEER® Light and SEER® Light Extend

Specifications available on request

SEER® Light Leadwire Set Kits – IEC			
Description	Length	Part Number	Quantity
Carry Pouch and 2 Channel Patient cable	140 cm/55 in	2042686-003	1
Carry Pouch and 3 Channel Patient cable		2042686-004	
Carry Pouch and 2 Channel Patient cable	60 cm/24 in	2042686-007	1
Carry Pouch and 3 Channel Patient cable		2042686-008	

Product compatibility

SEER® Light and SEER® Light Extend

Specifications available on request

Not available in all countries



SEER® Light Accessories		
Description	Part Number	Quantity
Compact Flash	2008750-003	1
Connect Kit	2020044-292	1
Controller Kit	2042685-007	1

Product compatibility

SEER® Light and SEER® Light Extend
Specifications available on request

SEER® Light Accessories		
Description	Part Number	Quantity
Connect Kit for China	2020044-301	1
Controller Kit for China	2042685-008	1

Sold ONLY in China

Product compatibility

SEER® Light and SEER® Light Extend
Specifications available on request



SEER® Light Holter Monitor Accessories		
Description	Part Number	Quantity
SEER® Light, shoulder straps, waist belt, portable pouch	2104827-001*	1

Product compatibility

2104827-001 is compatible with SEER® Light and SEER® Light Extend
Specification available on request

Not available in all countries





Product may not be available in all countries and regions.
Full product technical specifications are available upon request.
Contact a GE Healthcare representative for more information.

Please visit: www.gehealthcare.com/promotional-locations

Data subject to change.

GE Healthcare provides transformational medical technologies and services to meet the demand for increased access, enhanced quality and more affordable healthcare around the world. GE (NYSE: GE) works on things that matter – great people and technologies taking on tough challenges.

From medical imaging, software & IT, patient monitoring and diagnostics to drug discovery, biopharmaceutical manufacturing technologies and performance improvement solutions, GE Healthcare helps medical professionals deliver great healthcare to their patients.

© 2020 General Electric Company – All rights reserved.

GE, the GE monogram, Imagination at work, MAC™, CASE™, CardioSoft™, CardioMem™, SEER®, Aespire™, Aestiva™, Aisys™, Carestation™, CARESCAPE™, Solar™, Dash™, Entropy™, D-lite™, D-fend™, Corometrics™, Shuttle™, Panda™, Lullaby™ and TruSignal™ are trademarks of General Electric Company. ICU Medical Transpac is a registered trademark of ICU Medical, Inc.; Utah/Biotrans is a registered trademark of Utah Medical Products, Inc. Edwards Lifesciences TruWave is a trademark of Edwards Lifesciences Corporation. Argon is a trademark of Argon Medical Devices. AMSORB and G-Can are trademarks of Armstrong Medical. Hytrel is a registered trademark of E. I. du Pont de Nemours and Company or its affiliates. Inteliflo is a registered trademark of Flexicare Medical Ltd. Nellcor and Oximax are trademarks of a Medtronic Company. Masimo and SET are trademarks of Masimo Corporation. All other trademarks are the property of their respective owners.

Reproduction in any form is forbidden without prior written permission from GE. Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.

JB77272CA

Medical devices advertised in this catalogue may not be licensed in accordance with Canadian law.